

## Faecal Occult Blood Test

**Your Doctor has asked you to collect a sample of your faeces (poo) on 3 days for Pathology testing.**

Open the **brown topped container**.

Use the scoop attached to the lid to place a scoopfull of faeces in the container provided.

Close the lid on the container tightly.

Write your full name, date of birth, the date and time of the collection on each container.

Place container in specimen bag and seal.

Repeat this procedure once a day for 3 consecutive days.

Do not perform this collection if you have visible blood in your faeces or urine.

Store samples in fridge if samples cannot be delivered to the collection center daily.

Ensure samples and request slips are returned to Austin Pathology within 3 days of collection.

[www.austinpathology.org.au](http://www.austinpathology.org.au)

